

EC. STEAM 98

VACUUM TEST PACK BOWIE & DICK TYPE

- Disponible Bowie & Dick Type Autoclave Test Pack
- Paquet-test de Bowie & Dick à usage unique, pour autoclave
- Bowie & Dick - Typ Einmaltestpaket
- Pacco monouso per test di Bowie & Dick
- Paquete Prueba Descartable para Autoclave Tipo Bowie & Dick
- Bowie & Dick disposable autoclave-testpakket
- Bowie & Dick testpaket av engångstyp
- Bowie & Dick erangas testpakke
- Bowie & Dick erangas testpakke
- Kertakäyttöinen Bowie & Dick Autoklaavin testipakkaus
- Pacote prova descartável Bowie & Dick para autoclave

- Jednorazowy zestaw testowy typu Bowie & Dick
- Egyezkeshasználatos Bowie & Dick tesztsomag
- Jednorázový Bowie & Dick testovací balíček do autoclávu
- Jednorazowy Bowie & Dick test
- Testm paket tipa Bowie & Dick za enkratno uporabo
- Bowie & Dick testpakket
- Vieneiz letojams Bowie & Dick autoclavu paldaudes komplekts
- Vienekartinis Bowie & Dick tipo autoclavui testavimo paketas
- Tek kullatintlik Bowie & Dick Ottiklav test Paketi
- مستلزمك لاختبار القابلة بوفتي و نيك وحدة كاشية

AUTOCLAVE A VAPORE / STEAM AUTOCLAVE
AUTOCLAVE À VAPEUR / AUTOCLAVE DE VAPOR
C/C/D / CYCLE 134°C / 3.5 min



CE

Solo per
Only for

FR

TR

CZ



SISTEMI PER IL CONTROLLO DELLA STERILIZZAZIONE
Rif. NORM EN ISO 11140-1, CLASSE 2; EN ISO 11140-3 & 11140-4
ECS srl - ITALY - www.ecssrl.com - info@ecssrl.com